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
BIOQUANT Cortical Bone Report			
Animal ID:		Age:	
Animal Group:		Sex:	
Derived Indices		Value	Unit
Bright Field Low Power Measurements			
Ct B.Ar/T.Ar.	Cortical Bone Area to Tissue Area	67.67	%
Ec.Ar/T.Ar.	Endocortical Area to Tissue Area	32.33	%
	2nd Moment of Cortical Area	29.13	mm ⁴
Ps.Pm	Periosteal Perimeter	54.84	mm
Mean Ct.Wi.	Mean Cortical Width	2.91	mm
Anterior Ct.Wi.	Anterior Cortical Width	2.69	mm
Posterior Ct.Wi.	Posterior Cortical Width	4.48	mm
Medial Ct.Wi.	Medial Cortical Width	3.11	mm
Lateral Ct.Wi.	Lateral Cortical Width	3.05	mm
Bright Field High Power Measurements - Specialty Surfaces			
Pm/B.Ar.	Endocortical Perimeter to Bone Area	0.34	mm ⁻¹
Pm/T.Ar.	Endocortical Perimeter to Tissue Area	0.23	mm ⁻¹
O.Pm / Ec.Pm.	Osteoid Perimeter to Endocortical Perimeter	2.53	%
E.Pm / Ec.Pm.	Erosion Perimeter to Endocortical Perimeter	1.41	%
Q.Pm / Ec.Pm.	Quiescent Perimeter to Endocortical Perimeter	96.06	%
Ob.Pm / Ec.Pm.	Osteoblast Perimeter to Endocortical Perimeter	0.00	%
Oc.Pm / Ec.Pm.	Osteoclast Perimeter to Endocortical Perimeter	1.41	%
O.Pm(Ob+) / Ec.Pm.	Osteoid Perimeter with Osteoblasts to Endocortical Perimeter	0.00	%
E.Pm(Oc+) / Ec.Pm.	Erosion Perimeter with Osteoclasts to Endocortical Perimeter	1.41	%
Bright Field High Power Measurements - Osteon Characterization			
Mean On.Pm.	Mean Osteon Perimeter	416.71	mcm
Mean H.Pm.	Mean Haversian Perimeter	149.13	mcm
Mean Min.On.Dm.	Mean Minimum Osteon Diameter	98.79	mcm
Mean Max.On.Dm.	Mean Maximum Osteon Diameter	138.46	mcm

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Summary

This template summarizes the primary data and derived indices¹ found in the standard trabecular bone histomorphometry template of BIOQUANT OSTEON II. This is not a complete list.

Easily Add or Customize Indices

BIOQUANT data collection and report is highly adaptable and may be immediately expanded to include other primary and derived indices. The list below is a common beginning place that is expanded or simplified to suit the needs of the laboratory.

Calculated Indices

Please check the calculated indices you would like to be included in your template.

Volume Measurements

- TV Tissue Volume (mm²)
- BV Bone Volume (mm²)
- Md.V Mineralized Volume (mm²)
- OV Osteoid Volume (mm²)
- BV/TV Bone/Tissue Volume (%)
- Md.V/BV Mineralized/Bone Volume (%)
- Md.V/TV Mineralized/Tissue Volume (%)
- OV/BV Osteoid/Bone Volume (%)
- OV/TV Osteoid/Tissue Volume (%)
- Additional Unlimited Custom Indices

Thickness Measurements

- O.Th Osteoid Thickness (mcm)
- W.Th Wall Thickness (mcm)
- Ir.L.Th Interlabel Thickness (mcm)
- Additional Unlimited Custom Indices

Number Measurements

- N.Ob Osteoblast Number
- N.Oc Osteoclast Number
- N.Ob/OS Osteoblast Numb/Osteoid Surf (Ob/mm)
- N.Oc/ES Osteoclast Numb/Erosion Surf (Oc/mm)
- Additional Unlimited Custom Indices

Surface Measurements

- BS Bone Surface (mm)
- OS Osteoid Surface (mm)
- Ob.S Osteoblast Surface (mm)
- ES Erosion Surface (mm)
- Oc.S Osteoclast Surface (mm)
- QS Quiescent Surface (mm)
- sLS Single Labeled Surface (mm)
- dLS Double Labeled Surface (mm)
- BS/BV Bone Surface / BV (mm⁻¹)
- BS/TV Bone Surface / TV (mm⁻¹)
- OS/BS Osteoid/Bone Surface (%)
- ES/BS Erosion/Bone Surface (%)
- QS/BS Quiescent/Bone Surface (%)
- Ob.S/BS Osteoblast/Bone Surface (%)
- Oc.S/BS Osteoclast/Bone Surface (%)
- Additional Unlimited Custom Indices

Structural Indices

- Tb.Th Trabecular Thickness (Plate Model) (mcm)
- Tb.N Trabecular Number (Plate Model)(mcm⁻¹)
- Tb.Sp Trabecular Separation (Plate Model) (mcm)
- Additional Unlimited Custom Indices

Kinetic Indices

- MS Mineralizing Surface (mm)
- MS/BS Mineralizing/Bone Surface (%)
- MAR Mineral apposition rate (mcm/day)
- Aj.AR Adjusted apposition rate (mcm/day)
- BFR/BS Bone formation rate (mcm/day)
- Mlt Mineralization lag time (days)
- Omt Osteoid maturation time (days)
- FP Formation period (days)
- Rs.P Resorption period (days)
- Rv.P Reversal period (days)
- Rm.P Remodeling period (days)
- QP Quiescent period (days)
- Tt.P Total period (days)
- Ac.f Activation frequency (/year)
- Additional Unlimited Custom Indices

1. Parfitt AM, Drezner MK, Glorieux FH, Kanis JA, Malluche H, Meunier PJ, Ott SM, Recker RM. 1987. Bone histomorphometry: standardization of nomenclature, symbols, and units. J Bone Miner Res 2:595-610.

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Summary

This template summarizes the primary data and derived indices¹ found in the standard cortical bone histomorphometry template of BIOQUANT OSTEO II. This is not a complete list.

Easily Add or Customize Indices

BIOQUANT data collection and report is highly adaptable and may be immediately expanded to include other primary and derived indices. The list below is a common beginning place that is expanded or simplified to suit the needs of the laboratory.

Calculated Indices

Please check the calculated Indices you would like to be included in your template.

Area Measurements

- (Ct.B.Ar) Cortical Bone Area
- (Ec.Ar) Endocortical Area
- 2nd Moment of Ct. Area
- (T.Ar) Tissue Area
- (Ct.B.Ar/T.Ar) % Cortical Bone
- (Ec.Ar/T.Ar) % Endocortical Area
- (Sa.Ct.Ar) Sampled Cortical Area
- (On.Ar) Osteon Area
- On.Ar/ Sa.Ct.Ar % Osteon Area
- (H.Ar) Haversian Area
- H.Ar/ Sa.Ct.Ar % Haversian Area
- (On.Fg.Ar) Osteon Fragment Area
- On.Fg.Ar/ Sa.Ct.Ar % Fragment Area
- (Po.Ar) Pore Area
- (Po.Ar/Sa.Ct.Ar) Porosity
- (Oc.Ar) Osteocyte Area
- Unlimited Additional Indices

Perimeter Measurements

- (Ps.Pm) Periosteal Perimeter
- (Tt.On.Pm) Total Osteon Perimeter
- Mean On.Pm Mean Osteon Perimeter
- (H.Pm) Total Haversian Perimeter
- (H.Pm) Mean Mean Haversian Perimeter
- Ec.Pm Endocortical Perimeter
- Ec.O.Pm Osteoid Perimeter
- Ec.E.Pm Erosion Perimeter
- Ec.Q.Pm Quiescent Perimeter
- Ec.Ob.Pm Osteoblast Perimeter
- Ec.Oc.Pm Osteoclast Perimeter
- Ec.Pm/B.Ar Ec.Pm / Bone Area
- Ec.Pm/T.Ar Ec.Pm / Tissue Area
- Unlimited Additional Indices

Width Measurements

- Mean Ct.Wi Mean Cortical Width
- Anterior Ct.Wi Anterior Cortical Width
- Posterior Ct.Wi Posterior Cortical Width
- Medial Ct.Wi Medial Cortical Width
- Lateral Ct.Wi Lateral Cortical Width
- Min.On.Dm Mean Min. Osteon Diameter
- Max.On.Dm Mean Max. Osteon Diameter
- (W.Wi) Mean Wall Width
- Unlimited Additional Indices

Number Measurements

- (N.On) Osteon Number
- (OPD) Osteon Population Density
- (N.On.Fg) Osteon Fragment Number
- (Oc.N) Osteocyte Number
- (Oc.Dn) Osteocyte Density
- Unlimited Additional Indices

Fluorescent Indices

- (Ec.M.Pm) Ec. Mineralizing Pm
- (Ec.Ir.L.Wi) Avg Interlabel Wi
- (Ec.MAR) Mineral Apposition Rate
- (Ps.Ir.L.Wi) Avg Interlabel Wi
- (Ps.MAR) Mineral Apposition Ra
- (Ps.sL.Pm) Periosteal Single Label Pm
- (Ps.dL.Pm) Periosteal Double Label Pm
- (Ps.M.Pm) Periosteal Mineralizing Pm
- Unlimited Additional Indices

Formation Indices

- (Ac.F) Activation Frequency
- (BFR) Bone Formation Rate
- Unlimited Additional Indices

1. Parfitt AM, Drezner MK, Glorieux FH, Kanis JA, Malluche H, Meunier PJ, Ott SM, Recker RM. 1987. Bone histomorphometry: standardization of nomenclature, symbols, and units. J Bone Miner Res 2:595-610.

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